

# **B-Scan Plus/UBM Plus Recommended Computer Specifications Software Version 4.07.01**



## **Processor**

Intel Core 2 Duo 2.20 GHz, Core i3 2.2 GHz or higher

## **Memory**

3.0 GB, DDR2-667 SDRAM or higher

## **Video Card or Chip**

Desktop PC: Any NVIDIA or ATI Video Card with 256MB graphic memory or more.  
Laptop PC: Any NVIDIA, ATI or Intel Video Chip

## **Compatible Screen Resolutions**

### **1. For monitors 15.4 inches or larger: (recommended)**

Minimum supported resolution: 1366 X 768 pixels

Maximum supported resolution: 2560X 1440 pixels

*Use only the native resolution of your tablet or PC when using the B-scan/UBM Software.*

### **2. For tablets or monitors smaller than 15.4 inches:**

Minimum supported resolution: 1366 X 768 pixels

Maximum supported resolution: Software version 4.07.01 is compatible with previously unsupported higher resolutions. In order to enable proper functionality of these resolutions, your Windows display size setting may need to be adjusted. This setting can be accessed in Windows 7, 8, and 8.1 by navigating to “Control Panel”→ “Display” and selecting the correct display size setting for your native resolution. The display size setting compatible for your resolution is shown below:

<b>PC/Tablet Native Resolution</b>	<b>Windows 7, 8, 8.1 Display Size Setting</b>
Less than 1920 X 1080	100% (Default setting)
1920 X 1080	150%
1920 X 1200	150%
2160 X 1440	150%
2560 X 1440	150%

*Use only the native resolution of your tablet or PC when using the B-scan/UBM Software.*

# **B-Scan Plus/UBM Plus Recommended Computer Specifications Software Version 4.07.01**



## **Hard Disk**

Minimum 80 GB, 7200 RPM, recommended 160 GB or higher for high patient flow site (more than 30 B-Scans per month).

## **Optical Drive**

8X (or higher) DVD+/-RW with control software (needed for archiving)

## **USB Ports**

Minimum of 4 USB 2.0

## **Medical Grade Power Supply**

### **The following devices are not currently supported by the B-scan and UBM software:**

1. Lenovo Yoga Pro 2 Tablet (Has native resolution of 3200 X 1800)
2. Tablets or PCs with native resolutions above 2560 X 1440.
3. Tablets or PCs that do not use Microsoft Windows.

### **Software and Operation System (OS)**

Microsoft Word<sup>®</sup> 2003, 2007 and 2010 (for report generation)  
Windows 8<sup>®</sup> or 7<sup>®</sup>, Microsoft Windows<sup>®</sup> Vista or Microsoft Windows<sup>®</sup> XP Professional Edition, SP2, Microsoft Windows<sup>®</sup> 8.1

**If the computer on which you intend to load the B-Scan/UBM Software does not meet or exceed the listed specifications, we highly recommend that you do not install the B-Scan Plus/UBM Plus Software.**

**It is recommended that the computer used in conjunction with the B-Scan/UBM Plus system be dedicated for the sole use of Accutome ultrasound devices. Installation of software not vital to the operation of the B-Scan, UBM or other Accutome devices is highly discouraged, as it may detrimentally affect the performance of the B-Scan/UBM Plus system. Additionally, accessing the Internet may negatively affect the performance of the B-Scan/UBM and make the system vulnerable to viruses and other harmful external factors.**